


www.cowellmedi.com



DiaBONE

Bovine Bone Substitute



CWM

Cowellmedi Co., Ltd.


The pioneers in Dental Implant & E.rhBMP-2

DiaBONE


Bovine Bone Substitute

A Bone 100% fused to Natural Human Bone


- Fast blood penetration.
- Great hydrophilicity.
- 3-Dimensional structure.
- Easy to handle.
- Maximizes bone fusion.
- Mutually connected porous structure.
- Optimal cell attachment and blood absorption.
- Stimulates activity of osteoclast and osteoblast.



Human bone structure




Dia bone structure (x50)



Dia bone structure (x1,500)



Hydrophilicity 1




Hydrophilicity 2

Safe & Trustable Material


- Made of 100% bovine bone.
- Cleansing over 30 times to perfectly remove organic substances.
- Firmed bone formation as highly dense.
- 100% pure HA & 99.73% of bone crystallization.



Raw material



Graft test 1



Graft test 2

(New bone formation clearly observed around grafted bone site)

Volume and Particle Size

Product Code	Particle Size	Volume
G2015	0.25~1.0mm	0.15g
G2025	0.25~1.0mm	0.25g
G2050	0.25~1.0mm	0.5g
G2100	0.25~1.0mm	1g
G2200	0.25~1.0mm	2g

Product Code	Particle Size	Volume
G5015	1.0~2.0mm	0.15g
G5025	1.0~2.0mm	0.25g
G5050	1.0~2.0mm	0.5g
G5100	1.0~2.0mm	1g
G5200	1.0~2.0mm	2g

DiaBONE

Cowellmedi Co.,Ltd

Floor 6, Blue Fin Tower, 42, Seochojungang-ro, Seocho-gu, Seoul, Korea
Tel. +82-2-3453-5085 Fax. +82-2-3453-5086 E-mail. cib@cowellmedi.com

Cowellmedi USA INC

218 Trianon LN Villanova PA 19085-1442 USA
Tel. 1-623-939-1344 Fax. 1-623-939-1472

Cowell R&D Institute

48, Hakgam-daero 221beon-gil, Sasang-gu, Busan, 46986, Korea
Tel. +82-51-314-2028 Fax. +82-51-314-2026



www.cowellmedi.com